

V7 Portable Electric Suction Units



High quality low noise electric suction units for wards, theatre, recovery, ICU and A&E, plus a battery powered model for use in emergency vehicles

- ◆ Choice of power options: mains or mains plus battery
- ◆ Cost-effective solution for all clinical areas
- ◆ Collection jar fittings are compatible with other popular cannister/liner systems
- ◆ Integrated overflow valve
- ◆ Oil-free, no maintenance pump
- ◆ Bacteria filter and suction hose
- ◆ V7 M is supplied with a 12/24 v power cable for in-vehicle use

Partnership for Life

Additional Features

- ◆ Robust ergonomic case with integrated handle
- ◆ Oil-free, no-maintenance pump
- ◆ Manometer scale from 0 – 100 kPA (0 – 760 mmHg)
- ◆ Vacuum control knob
- ◆ 84 kPa (630 mmHg) maximum vacuum
- ◆ 30 L/min suction capacity
- ◆ Low noise level: 49 ±1.5 dB (AS) 1 m
- ◆ Bacteria filter and suction hose included
- ◆ Water resistant ON/OFF switch. Battery models are fitted with an additional switch protector and a LED display (V7 M)
- ◆ Fittings compatible with other popular canister and disposable liner systems

Technical Specification

	V7 Plus AC	V7 M
Weight	3.5 kg	5.2 kg
Power Supply	220 VAC	110/220 VAC, 12/24 V (DC Battery)
Power Consumption	120 VA	130 VA (AC), 84 W (DC)
Rechargeable Battery	-	12 V (4.5 Ah) Lead Acid
Use with Battery	-	65 minutes
AC Internal Battery Charger	-	Yes
12/24 V Cable (for Ambulance use)	-	Yes
Bacteria Filter	Yes	Yes
Suction Capacity	30 L/min	30 L/min
Maximum Vacuum	-84 kPa (-630 mmHg)	-84 kPa (-630 mmHg)

Ordering Information and Accessories

V7 Plus AC

Part No.	Description
54050	V7 Plus AC Suction Unit
54051	V7 Plus AC Suction Unit with one litre collection jar

V7 M

Part No.	Description
54053	V7 M Suction Unit
54054	V7 M Suction Unit with one litre collection jar

Accessories

Part No.	Description
54056	Rolling stand
54057	Transport bag
54058	Bacteria filter

Suction Technology

Distributed by: Penlon Limited • Manufactured by: Hersill, S.L.



Penlon Limited, Abingdon Science Park, Barton Lane, Abingdon, OX14 3PH UK

UK

Tel: 01235 547036
 Fax: 01235 547023
 E-mail: uksales@penlon.co.uk

International

Tel: +44 1235 547001
 Fax: +44 1235 547021
 E-mail: export@penlon.co.uk



Penlon Limited reserves the right to change technical specification without prior notice.

© 2009 Penlon Limited

Penlon Limited is a Member of the InterMed Group