


Penlon SP M8 Patient Monitor

ANAESTHESIA SOLUTIONS

- ◆ Flexible solutions to suit changing clinical situations
- ◆ Five on-board slots and optional module rack with eight additional slots
- ◆ Up to twelve waveforms
- ◆ Easy to upgrade
- ◆ 17" colour display with optional touch screen



Penlon SP M8 Patient Monitor



Features and options configured to your exact specification

- 1 17" colour display with optional touch screen
- 2 Integrated carry handle
- 3 High visibility, prioritized alarms
- 4 Multi-parameter module (MPM)
- 5 Intelligent, quiet cooling system
- 6 Interface control prompt
- 7 Configurable QuickKeys for fast access to functions
- 8 Up to eight simultaneous waveforms

The Penlon SP M8 patient monitor provides clarity, simplicity, flexibility and fast processing power for superior monitoring in all clinical areas.

Ultimate flexibility begins with multiple module options for a wide variety of clinical situations, including the operating theatre and recovery.

Parameters and modules

ECG

- 3/5 lead or optional 12 lead set
- ST detection and arrhythmia analysis
- Optional 'extended' Mortara algorithms

SpO₂

Choice of SpO₂ including optional Masimo SET™ or Nellcor OxiMax™.

IBP

Up to eight IBP channels.

Gases

Mainstream, Sidestream or Microstream CO₂ modules, Multi-gas with optional paramagnetic O₂ and integrated Bispectral Index (BIS).

Additional parameters

- Cardiac output

Multi-parameter module

The multi-parameter module (MPM) measures SpO₂, NIBP, 3/5 lead ECG, Temperature and Respiration as standard. There are eleven further MPM configurations to choose from.

Multi-gas for anaesthesia

The SP M8 is easily configured for use in the operating theatre with a choice of automatic gas identification modules, ensuring continuous monitoring of O₂, CO₂, N₂O, and anaesthetic agent. Depth of anaesthesia can be evaluated by the integral BIS unit.

Display

The 17" touch screen display is fast and easy to use, and can be configured to show up to eight simultaneous waveforms, with user-configured waveform colours.

1 Minitrends display

Split screen view showing mintrends for multiple parameters alongside waveforms.

2 Large numerics display

Shows vital information clearly, to minimise interpretation errors.

3 OxyCRG display

Combines multiple respiratory trends in a single area, allowing instant evaluation of breathing efficiency.

4 Menu

QuickKeys give instant access to basic tasks, including patient admission and discharge, and alarm configuration.


Networking and data storage

Wired or wireless networking for connection to central station care groups (up to 10 monitors per group) allows viewing of other beds. Additional local data storage is available via the optional CF card.

Maintenance and after-sales support

Penlon is committed to a successful, long term relationship with all our customers. Comprehensive warranty provides user peace of mind and after-sales support.

Additional services and warranties can be purchased to meet your particular needs.




About Penlon ♦

Penlon was founded in 1943 by personnel from the Department of Anaesthesia at Oxford University. One of the first products was the Macintosh Laryngoscope, then a revolutionary design, and still the most widely used today, invented by the late Sir Robert Macintosh, Professor of Anaesthetics at Oxford University.

Today Penlon continues to design, engineer and build high quality anaesthesia products at its UK headquarters. The company is proud to have over 70 years' dedicated experience, many awards for product design, and an impressive four Queen's Awards for Enterprise, one for 'Innovation' and three for 'International Trade'.

Penlon devices feature intuitive user interfaces that require minimal operator training, putting clinicians in control, enabling them to focus on what is most important – patient safety and wellbeing.



 Global Headquarters  Local Office  Distribution Partner



Like us at facebook.com/PenlonGlobal



Follow us at twitter.com/PenlonGlobal

0088

Doc. No. SPM8 0177SL / April 2017. Penlon and InterMed are trademarks of Penlon Limited. All other trademarks are the property of their respective owners.

Part of the InterMed Group © Penlon Ltd. 2017 All rights reserved

Penlon Limited
Abingdon Science Park
Barton Lane, Abingdon
OX14 3NB, UK

General
t +44 (0) 1235 547000
f +44 (0) 1235 547041
w www.penlon.com

International Sales
t +44 (0) 1235 547001
f +44 (0) 1235 547021
e international.sales@penlon.com

UK Sales
t +44 (0) 1235 547036
f +44 (0) 1235 547023
e uk.sales@penlon.com

Technical Support
t +44 (0) 1235 547060
f +44 (0) 1235 547061
e tech.support@penlon.com

