

Penlon Sigma Delta MRI Vaporizer

ANAESTHESIA SOLUTIONS

- ◆ Approved for use with 1.5 and 3 tesla scanners, up to 1000 gauss¹
- ◆ Multiple agent and filler block options
- ◆ Superb performance, particularly at low flow
- ◆ Service Free²
- ◆ Low body weight



TESTED AND APPROVED FOR USE WITH
**GE Healthcare, Philips and
Siemens MRI Scanners³**

Penlon Sigma Delta MRI Vaporizer

ANAESTHESIA SOLUTIONS




The Sigma Delta MRI vaporizer delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flow


Vaporizer Type	Sevoflurane 8%			Isoflurane 5%		Halothane 5%
	Filler Type			Filler Type		Filler Type
	Pour fill	Key fill	Quik-Fit	Pour fill	Key fill	Key fill
Selectatec Compatible with Interlock	✓	✓	✓	✓	✓	✓
Cagemount Taper 23 mm	✓	✓	✓	✓	✓	✗
Dräger Auto-Exclusion® Compatible ⁴	✓	✓	✓	✓	✓	✗

- ◆ Key fill, Quik Fil[®], or Pour fill
- ◆ Sevoflurane, Isoflurane and Halothane
- ◆ Selectatec[®], Cagemount and Dräger Auto-Exclusion^{®4}
- ◆ Internationally recognised, award-winning design
- ◆ Clearly labelled as MR Conditional


Typical performance graph




Backbar mounting systems



Selectatec Compatible with Interlock



Cagemount Taper 23 mm



Dräger Auto-Exclusion Compatible⁴

CE 1639 UK CA 0120



Doc. No. DM0123SL / November 2023. Penlon and Sigma Delta are trademarks of Penlon Limited. All other trademarks are the property of their respective owners. ¹ Conditions apply. ² Presumes ten year product life requiring no preventive maintenance service. It is recommended that a service is carried out at ten years. Halothane vaporizers require a five year service. ³ Visit www.penlon.com for further information and to download compatibility statements. ⁴ Auto-Exclusion not available in the US.

Rx only



A Medcaptain Company

© Penlon Ltd. 2023 All rights reserved

Penlon Limited
Abingdon Science Park
Barton Lane, Abingdon
OX14 3NB, UK

General
t +44 (0) 1235 547000
w www.penlon.com

International Sales
t +44 (0) 1235 547001
e international.sales@penlon.com

UK Sales
t +44 (0) 1235 547036
e uk.sales@penlon.com

Technical Support
t +44 (0) 1235 547060
e tech.support@penlon.com

EC REP

OBELIS S.A.
Bd Général Wahis, 53
B-1030 Brussels
Belgium