

Penlon Sigma Delta Vaporizer

ANAESTHESIA SOLUTIONS

- ◆ Multiple agent and filler block options
- ◆ Five backbar mounting systems, including Dräger Auto-Exclusion®†
- ◆ Superb performance, particularly at low flow
- ◆ Service Free*
- ◆ Low body weight



Penlon Sigma Delta Vaporizer

ANAESTHESIA SOLUTIONS




The Sigma Delta vaporizer delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flow


Vaporizer Type	Sevoflurane 8%			Isoflurane 5%		Halothane 5%	
	Filter Type			Filter Type		Filter Type	
	Pour fill	Key fill	Quik-Fit	Pour fill	Key fill	Pour fill	Key fill
Selectatec Compatible with Interlock	✓	✓	✓	✓	✓	✓	✓
Dräger Plug-In Interlock Compatible	✓	✓	✓	✓	✓	✓	✓
Dräger Auto-Exclusion® Compatible*	✓	✓	✓	✓	✓	✓	✓
Cagemount Taper 23 mm	✓	✓	✓	✓	✓	✓	✓
North American Dräger Compatible	✓	✓	✓	✓	✓	✓	✓

- ◆ Key fill, Quik Fil®, or Pour fill
- ◆ Sevoflurane, Isoflurane and Halothane
- ◆ Selectatec®, Dräger Plug-In®, Dräger Auto-Exclusion®+, North American Dräger, Cagemount
- ◆ Internationally recognised, award-winning design

Typical performance graph



Backbar mounting systems



facebook.com/PenlonGlobal



linkedin.com/company/penlon-ltd



twitter.com/PenlonGlobal



1639 Doc. No. SD0721SL / October 2021. Penlon, InterMed and Sigma Delta are trademarks of Penlon Limited. All other trademarks are the property of their respective owners. *Presumes ten year product life requiring no preventive maintenance service. It is recommended that a service is carried out at ten years. Halothane vaporizers require a five year service.†Auto-Exclusion not available in the US.

Part of the InterMed Group © Penlon Ltd. 2021 All rights reserved

Penlon Limited
Abingdon Science Park
Barton Lane, Abingdon
OX14 3NB, UK

General
t +44 (0) 1235 547000
f +44 (0) 1235 547041
w www.penlon.com

International Sales
t +44 (0) 1235 547001
f +44 (0) 1235 547021
e international.sales@penlon.com

UK Sales
t +44 (0) 1235 547036
f +44 (0) 1235 547023
e uk.sales@penlon.com

Technical Support
t +44 (0) 1235 547060
f +44 (0) 1235 547061
e tech.support@penlon.com

EC REP

OBELIS S.A.
Bd Général Wahis, 53
B-1030 Brussels
Belgium