



A UK COMPANY



Accessories and Consumables

PATIENT MONITORING SOLUTIONS

Penlon 

The AnyView Patient Monitor accessory range is precisely engineered to meet the highest standards of quality, performance and reliability.



The extensive range including ECG, EEG and SpO₂ accessories, is the ideal choice for use across the hospital.

Contents





- 3 Standard and Optional Accessories
- 4 ECG
- 5 BLT SpO₂
- 6 NIBP
- 7 TEMP, IBP and Masimo Mainstream EtCO₂
- 8 C.O., BLT CO₂ and Nellcor SpO₂
- 9 Masimo SpO₂ and BLT Sidestream CO₂
- 10 Suntech NIBP and BIS
- 11 EEG and RM
- 12 Masimo Sidestream and Multi-Gas Mainstream



Standard Accessories



Product	Name	Part Number	Length	Description	Compatibility
	IEC integrated 3-lead snap	5010282	3.2 m	Adult/paediatric use, defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	Integrated SpO ₂ sensor	5010426	2.9 m	Adult, analog, soft reusable, encrypted	P: 1/12/15/18/22 S: 10/12 M: 10/12
NO IMAGE AVAILABLE	Power cord	5010277	-	EU standard	P: 12/15/18/22 S: 10/12
NO IMAGE AVAILABLE	Lithium battery	5010276	-	11.1 V / 2500 mAh	S: 10/12 M: 10/12
NO IMAGE AVAILABLE	Manual	5010278	-	User Manual	P: 12/15/18/22
NO IMAGE AVAILABLE	5-lead snap	5010284	-	Connector ECG cable (IEC) (12 pin) (AQS) (encrypted)	P: 12/15/18/22 S: 10/12
NO IMAGE AVAILABLE	Rechargeable lithium battery	5010374	-	11.1 V / 5000 mAh (P series)	P: 12/15/18/22

Optional Accessories

Product	Name	Part Number	Adaptor	Description	Compatibility
	Auxillary module slot	5010419	-	NC/AMS connectors are required Supports up to 8 modules	P: 15/18/22
NO IMAGE AVAILABLE	Lithium battery	5005260	-	11.1 V / 2500 mAh	S: 10/12 M: 10/12
	Premium trolley	5010423	-	Universal across compatible BLT devices (Pictures are for illustration purposes only, actual product may vary)	P1:(with DS)/P12 S: 10/12 M: 10/12
	Premium trolley	5010345	Quick lock	Universal across compatible BLT devices (Pictures are for illustration purposes only, actual product may vary)	P1:(with DS)/P12 S: 10/12 M: 10/12
	Wall mount	5010346	Quick lock	Connects P12/S10/S12/M10/M12 to non-quick lock mountings	P: 1/12 S: 10/12 M: 10/12

ECG








Product	Name	Part Number	Length	Description	Compatibility
	IEC integrated 12-lead snap	5010279	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	IEC Integrated 6-Lead Snap	5010281	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 S: 10/12 M: 10/12
	IEC integrated 5-lead snap	5010282	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	AHA integrated 12-lead snap	5010283	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	AHA integrated 5-lead snap	5010284	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	IEC integrated 3-lead snap	5010285	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	AHA integrated 3-lead snap	5010286	3.2 m	Adult/paediatric use, Defibrillator proof, 1 pc/package	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	AHA split 5-lead clamp	5010287	3.9 m	Anti-static, anti-cautery, includes trunk cable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	IEC split 5-lead clamp	5010288	3.9 m	Anti-static, anti-cautery, includes trunk cable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	3-lead trunk cable	5010289	2.5 m	For 3 leads wires only	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	AHA plug-in 3-lead snap	5010290	0.6 m	Neonatal use, Defibrillator proof, 1 pc/package Requires 3-lead trunk cable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	IEC plug-in 3-lead snap	5010291	0.6 m	Neonatal use, Defibrillator proof, 1 pc/package. *Requires 3-lead trunk cable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
NO IMAGE AVAILABLE	3-lead snap connector ECG cable	5010292	-	(AHA) [12 pin] (plug-in) [AQS]	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12

Product	Name	Part Number	Length	Description	Compatibility
	Adult	5010280	10 pcs	Room temperature storage, Pre-applied electrogel, 10 pc/pack	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 M: 800/850
	Paediatric, neonate	5010293	10 pcs	Room temperature storage, pre-applied electrogel, 10 pc/pack	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 M: 800/850



BLT SpO₂

Product	Name	Part Number	Length	Description	Compatibility
	Integrated SpO ₂ sensor	5010295	2.9 m	Adult, analog, clip Reusable, encrypted (12 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Integrated SpO ₂ sensor	Integrated SpO ₂ sensor	1.5 m	Adult, analog, soft Reusable, encrypted (12 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	SpO ₂ analog extension cable	5010296	2.2 m	12 pin to 9 pin Analog SpO ₂ sensor use only	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Plug-in SpO ₂ sensor	5010297	1.5 m	Adult, analog, clip type Reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12M 800/850/860/880
	Plug-in SpO ₂ sensor	5010298	1.5 m	Adult, analog, soft Reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12M
	Plug-in SpO ₂ sensor	5010299	1.5 m	Paediatric, analog, clip Reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Plug-in SpO ₂ sensor	5010300	1.5 m	Paediatric, Analog, soft Reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Plug-in SpO ₂ sensor	5010301	1.5 m	Neonate, analog, strap Reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Plug-in SpO ₂ sensor	5010302	1.5 m	Adult/neonate, adhesive analog, disposable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12
	Multi-site Y sensor	5010303	2.9 m	Adult/paediatric/neonate Analog, reusable, encrypted Requires SpO ₂ analog extension cable (9 pin)	P: 1/12/15/18/22 S: 10/12 M: 10/12






NIBP

Product	Name	Part Number	Length	Description	Compatibility
	NIBP Hose	5010304	2.5 m	PVC material	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6
	Adult Cuff	5010427	-	Adult, 25-35 cm PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9 M: 860
	Adult Cuff	5010305	-	Adult, Large, 33-47 cm PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9 M: 860
	Adult Cuff	5010306	-	Adult, Thigh, 44-53 cm PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9
	Pediatric Cuff	5010307	-	Paediatric, 18-26 cm PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9 M: 860
	Neonatal Cuff	5010310	-	Neonatal, 6-11 cm PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9 M: 860
	Neonatal cuff	5010311	-	Neonatal, disposable PU material, latex free	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6/9
	Neonatal cuff	5010312	-	Neonatal, 3-5.5 cm Disposable PU material, latex free	
	Neonatal cuff	5010313	-	Neonatal, 4-8 cm Disposable PU material, latex free	
	Neonatal cuff	5010314	-	Neonatal, 6-11 cm Disposable PU material, latex free	
	Neonatal cuff	5010315	-	Neonatal, 7-13 cm Disposable PU material, latex free	


TEMP

Product	Name	Part Number	Length	Description	Compatibility
	Skin probe	5005296	3 m	Reusable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12
	Rectal probe	5005295	3 m	Reusable	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12





IBP

Product	Name	Part Number	Length	Description	Compatibility
	IBP extension cable	5010319	0.7 m	4 pin to 6 pin Reusable A/Q series only	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	Utah transducer cable	5010320	3.2 m	Require IBP extension cable	P: 1/12/15/18/22 Q: 3/5/7 S: 12 Utah transducer
	Utah IBP transducer set	5010321	-	IBP extension cable and Utah extension cable required	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	BD transducer cable	5010322	3.2 m	Require IBP extension cable	P: 1/12/15/18/22 Q: 3/5/7 S: 12 BD transducer
	BD IBP transducer set	5010323	-	IBP extension cable and BD extension cable required	P: 1/12/15/18/22 Q: 3/5/7 S: 12





Masimo Mainstream EtCO₂

Product	Name	Part Number	Length	Description	Compatibility
	Masimo gas module extension cable	5010332	2.2 m	Connect to: Masimo IRMA EtCO ₂ sensor, Masimo ISA EtCO ₂ sensor (for S12), Masimo IRMA AX+ Multi-Gas sensor	P: 1/12/15/18/22 Q: 3/5/7 S: 12


C.O.

Product	Name	Part Number	Length	Description	Compatibility
	Cardiac output cable	5010326	3 m	6 pin Edwards interface	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	In-line injection temperature probe	5010327	-	BD Connects to 5010328 for measuring temp of injectant	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	In-line injection temperature probe housing	5010328	3 m	BD Connects to 5010327 and catheter	P 1/12/15/18/22 Q 3/5/7 S 12
	Control syringe	5010329	-	12 ml Connects to 5010328	P: 1/12/15/18/22 Q: 3/5/7 S: 12


BLT CO₂

Product	Name	Part Number	Length	Description	Compatibility
	EtCO ₂ module extension cable	5010332	0.6 m	BLT CO ₂ extension cable required 50 ml/min sampling rate	S: 12 V: 6/9
	BLT T4F water filter	5010339	-	Disposable BLT Microflow CO ₂ module required	P: 1/12/15/18/22 S: 12 V: 6/9 M: 880
	Sampling line	5010333	2.4 m	Disposable CO ₂ sampling line	P: 1[DS required] P: 1/12/15/18/22 Q: 3/5/7 S: 12 V: 6/9 M: 880
	BLT CO ₂ nasal cannula	5010341	-	Adult/paediatric Disposable CO ₂ nasal cannula	P: 1/12/15/18/22 S: 12 V: 6/9 M: 880



Nellcor SpO₂

Product	Name	Part Number	Length	Description	Compatibility
	Reusable Clip	5010388	1 m	Adult, weight >40 kg Nellcor SpO ₂ adapter cable required	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9 M: 800/850/860/880







Masimo® SpO₂

Product	Name	Part Number	Length	Description	Compatibility
	RD Reusable sensor clip	5010375	3 ft	Adult, Weight >30 kg, Low perfusion, anti-motion 1 pc/box Requires RD extension cable	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9
	RD Reusable sensor clip	5010376	3 ft	Paediatric, 10-50 kg, Analog, 1 pc/box Requires RD extension cable	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9
	RD Adhesive sensor	5010379	-	Neonate (<3 kg) or Adult (>40 kg), Analog, adhesive wrap, 20 pcs/box Require RD extension cable	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9
	RD Extension cable	5010380	12 ft	Compatible with Masimo and Masimo Rainbow SET	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9
	Silicone wraps	5010381	-	Reusable silicone wraps for 5010357	P 1/12/15/18/22 Q 3/5/7 V 6/9
	RD Multi-site sensor	5010357	3 ft	Adult/Pediatric/Neonate Weight >1 kg, Analog, Multiple foam and adhesive wrap, 1 pc/box Require RD extension cable	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9
	RD Adhesive sensor	5010378	-	Infant, weight 3-20 kg, Analog, adhesive wrap, 20 pcs/box Require RD extension cable	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9




BLT Sidestream CO₂

Product	Name	Part Number	Length	Description	Compatibility
	BLT sidestream water trap	5010411	-	Disposable	P: 1(DS required) P: 12/15/18/22 Q: 3/5/7
	BLT sampling line	5010333	2.4 m	Disposable	P: 1(DS required) P: 12/15/18/22 Q: 3/5/7 S: 12 V: 6/9 M: 880
NO IMAGE AVAILABLE	-	5010340	-	L type 3 way airway adapter for intubated patient	-





Suntech NIBP

Product	Name	Part Number	Length	Description	Compatibility
	NIBP hose	5010304	2.5 m	PVC material	P: 1/12/15/18/22 Q: 3/5/7 S: 10/12 M: 10/12 V: 6
NO IMAGE AVAILABLE	Adult cuff	5010392	1 m	Adult cuff with connector, 23-33 cm	P: 1/12/15/18/22 Q: 3/5/7 V: 6/9 M: 800/850/860/880
	Neonatal cuff	5010393	3-6 cm	Neonatal, disposable	P: 1/12/15/18/22 Q: 3/5/7
	Neonatal cuff	5010394	4-8 cm	Neonatal, disposable	P: 1/12/15/18/22 Q: 3/5/7
	Neonatal cuff	5010395	6-11 cm	Neonatal, disposable	P: 1/12/15/18/22 Q: 3/5/7
	Neonatal cuff	5010396	7-13 cm	Neonatal, disposable	P: 1/12/15/18/22 Q: 3/5/7
	Neonatal cuff	5010397	8-15 cm	Neonatal, disposable	P: 1/12/15/18/22 Q: 3/5/7



BIS

Product	Name	Part Number	Length	Description	Compatibility
	BIS adult quarto sensor	5010399	-	Disposable, 1 set/packet	P: 1/12/15/18/22 Q: 3/5/7
	BIS adult extended use sensor	5010400	-	Disposable, 1 set/packet ICU use (Measures brain activity in adult patients requiring longer periods of monitoring)	P: 1/12/15/18/22 Q: 3/5/7
	BIS paediatric sensor	5010401	-	Disposable, 1 set/packet	P: 1/12/15/18/22 Q: 3/5/7




EEG

Product	Name	Part Number	Length	Description	Compatibility
NO IMAGE AVAILABLE	EEG accessories kit	5010402	-	EEG accessories kit (EEG patient cable 1 set, EEG cap 3 pcs/set, EEG Coupling agent 1 pc, syringe 1 set)	-
	EEG patient	5010403	-	For adult/paediatric use	P: 12/15/18/22 Q: 3/5/7 (Non-EU)
	EEG Cap	5010404	-	Adult/paediatric Small, 50c m	P: 12/15/18/22 Q: 3/5/7 (Non-EU)
	EEG Cap	5010405	-	Adult/paediatric Medium, 55 cm	P: 12/15/18/22 Q: 3/5/7 (Non-EU)
	EEG Cap	5010406	-	Adult/paediatric Large, 60 cm	P: 12/15/18/22 Q: 3/5/7 (Non-EU)
NO IMAGE AVAILABLE	EEG GEL	5010407	-	EEG Coupling Agent (473G)	-
NO IMAGE AVAILABLE	EEG Syringe	5010408	-	EEG Coupling Agent Syringe 1 set	-


RM

Product	Name	Part Number	Length	Description	Compatibility
	Flow sensor	5010409	3.3 m	Adult/paediatric	P: 12/15/18/22
	Flow sensor	5010410	3.3 m	Neonate	P: 1/12/15/18/22 Q: 3/5/7

Masimo® Sidestream

Product	Name	Part Number	Length	Description	Compatibility
	Masimo Nomoline HH nasal cannula	5010413	2 m	For non-intubated patients For adult patients For long-term applications in high humidity conditions	P: 12/15/18/22 Q: 3/5/7
	Masimo Nomoline HH airway adapter set	5010414	3 m	For adult and paediatric patients For long-term applications in high humidity conditions	P: 12/15/18/22 Q: 3/5/7
	Masimo NomoLine HH infant nasal cannula	5010415	2 m	For non-intubated patients For long-term applications in high humidity conditions	P: 12/15/18/22 Q: 3/5/7
NO IMAGE AVAILABLE	Masimo Nomoline HH airway adapter set	5010416	-	NomoLine HH (High Humidity) Airway Adapter Set Infant. (3829)	P: 12/15/18/22 Q: 3/5/7

Multi-Gas Mainstream


Product	Name	Part Number	Length	Description	Compatibility
	IRMA AX+ Multi-Gas sensor	5010418	-	IRMA AX+ Multi-Gas sensor (Main Stream)	-
	Airway adapter - neonatal/infant	5010335	-	Disposable airway adapter - Neonatal/Infant (4500)	-
	Airway adapter - adult/paediatric	5010336	-	Disposable airway adapter - Adult / paediatric (4499)	P: 1/12/15/18/22 Q: 3/5/7 S: 12
	BLT CO ₂ extension cable	5010332	2.2 m	Connects to: BLT Mainstream EtCO ₂ sensor BLT Microflow EtCO ₂ sensor	S: 12 V: 6/9

About Penlon ♦

Penlon was founded in 1943 by personnel from the Department of Anaesthesia at Oxford University. One of the first products was the Macintosh Laryngoscope, then a revolutionary design, and still the most widely used today, invented by the late Sir Robert Macintosh, Professor of Anaesthetics at Oxford University.

Today Penlon continues to design, engineer and build high quality anaesthesia products at its UK headquarters. The company is proud to have 80 years' dedicated experience, many awards for product design, and an impressive four Queen's Awards for Enterprise, one for 'Innovation' and three for 'International Trade'.

Penlon devices feature intuitive user interfaces that require minimal operator training, putting clinicians in control, enabling them to focus on what is most important – patient safety and wellbeing.



PenlonGlobal



company/penlon-ltd



PenlonGlobal



PenlonGlobal

Doc. No. AC0125SL Rev. A / 17 March 2025. Penlon is a trademark of Penlon Limited. All other trademarks are the property of their respective owners.



A Medcaptain Company

© Penlon Ltd. 2025 All rights reserved. Distributed by: Penlon Limited • Manufactured by: Biolight Co., Ltd.

Penlon Limited
Abingdon Science Park
Barton Lane, Abingdon
OX14 3NB, UK

General
t +44 (0) 1235 547000
w www.penlon.com

International Sales
t +44 (0) 1235 547001
e international.sales@penlon.com

UK Sales
t +44 (0) 1235 547036
e uk.sales@penlon.com

Technical Support
t +44 (0) 1235 547060
e tech.support@penlon.com

EC REP

Shanghai International
Holding Corp. GmbH (Europe)
Eiffestraße 80,
20537 Hamburg Germany