

# Penlon MRI Green System Laryngoscope Range

## INTUBATION SOLUTIONS

- ◆ Designed for use in the MR environment
- ◆ Safe at the 1000 gauss line
- ◆ Tested in 1.5 and 3.0 Tesla MRI scanners
- ◆ Macintosh, Miller, Robertshaw, Seward and McCoy blade styles
- ◆ Standard size handle
- ◆ High intensity LED light bulb
- ◆ MR-conditional batteries



# Penlon MRI Green System Laryngoscope Range

INTUBATION SOLUTIONS



## Use your Penlon MRI Green System laryngoscope safely in MR conditions with a static magnetic field of 1000 gauss/cm or less

The Penlon MRI Green System Laryngoscope Range has been developed through research into the requirements of both diagnostic and interventional Magnetic Resonance (MR) intubation.

Clinicians can intubate in the MR environment safe in the knowledge that the MRI Green System Range has been independently tested and certified to be MR-Conditional in a static magnetic field of 1000 gauss/cm or less.

### LED lamp

The laryngoscope's high-intensity LED lamp provides a clear, bright and cool light output via the fibre optic light guide, and is powered by replaceable MR-conditional alkaline batteries (C size).

### Non-magnetic battery

The MR-conditional alkaline batteries provides equivalent performance compared to standard batteries normally used in intubation instruments, and contains no lithium.

### Blades and handles

The fibre optic laryngoscope blades are available in a range of sizes in Macintosh, Miller, Robertshaw and Seward patterns, including the McCoy range. The laryngoscope handle is available in standard configuration.

The range of handle and blades means that clinicians in the MR environment have easy access to laryngoscopes suited to their requirements for the patient profile.

- ◆ Independently tested and certified to be MR-conditional in 1.5 and 3 Tesla MRI scanners
- ◆ Blade and handle combinations for the majority of patient profile groups
- ◆ LED lamp delivers a high intensity, cool, light beam
- ◆ MR-conditional alkaline batteries

- 1 Macintosh blades**  
Available in sizes 0, 1, 2, 3, 4 and 5
- 2 Miller blades**  
Available in sizes 00, 0, 1, 2, and 3
- 3 Robertshaw blades**  
Available in sizes 0 and 1
- 4 Seward blades**  
Available in sizes 1 and 2



[facebook.com/PenlonGlobal](https://www.facebook.com/PenlonGlobal)



[linkedin.com/company/penlon-ltd](https://www.linkedin.com/company/penlon-ltd)



[twitter.com/PenlonGlobal](https://www.twitter.com/PenlonGlobal)



Doc. No. MRG0121SL / June 2021. Penlon and InterMed are trademarks of Penlon Limited. All other trademarks are the property of their respective owners.

Part of the InterMed Group © Penlon Ltd. 2021 All rights reserved

Penlon Limited  
Abingdon Science Park  
Barton Lane, Abingdon  
OX14 3NB, UK

**General**  
t +44 (0) 1235 547000  
f +44 (0) 1235 547041  
w [www.penlon.com](http://www.penlon.com)

**International Sales**  
t +44 (0) 1235 547001  
f +44 (0) 1235 547021  
e [international.sales@penlon.com](mailto:international.sales@penlon.com)

**UK Sales**  
t +44 (0) 1235 547036  
f +44 (0) 1235 547023  
e [uk.sales@penlon.com](mailto:uk.sales@penlon.com)

**Technical Support**  
t +44 (0) 1235 547060  
f +44 (0) 1235 547061  
e [tech.support@penlon.com](mailto:tech.support@penlon.com)

EC REP

OBELIS S.A.  
Bd Général Wahis, 53  
B-1030 Brussels  
Belgium