

# Penlon Sigma EVA Suprane® (desflurane) Vaporizer

## ANAESTHESIA SOLUTIONS

- ◆ Lightweight, innovative, low life costs
- ◆ UK design and manufacture
- ◆ Patented technology to ensure optimum drug delivery
- ◆ Fast warm-up
- ◆ Compatible with Selectatec®, Dräger Plug-In® and Dräger Auto-Exclusion® mounting systems
- ◆ SAFE-FIL® filling system and high capacity reservoir



# Penlon Sigma EVA Suprane® (desflurane) Vaporizer

ANAESTHESIA SOLUTIONS



Accurate, lightweight, fast warm-up,  
low power consumption, low life costs.  
Maximum functionality, minimum  
expenditure.

- ◆ Approved for use with Suprane® (desflurane) by Baxter
- ◆ Intuitive design
- ◆ Carrier gas compensation
- ◆ Accurate agent delivery
- ◆ Low power consumption
- ◆ Up to 30 minutes battery back up
- ◆ Low agent indicator
- ◆ Prioritised alarms
- ◆ 1% – 18% concentration range
- ◆ 24 month warranty

## Innovative Desflurane Delivery Solution

The unique patented carrier gas compensation system ensures optimum drug consumption, delivering accurate agent output even at low flow rates. The vaporizer continuously compensates for changes in the O<sub>2</sub>, Air and N<sub>2</sub>O mixture set by the anaesthetist.

## Design and build

The lightweight, intuitive design of the Sigma EVA makes it easy to transport and use. The vaporizer can be fitted to the majority of anaesthesia machines in use today, and is compatible with Selectatec®, Dräger Plug-In®, and Dräger Auto-Exclusion® backbar mounting systems. The vaporizer is approved for use with Suprane® (desflurane) and can be quickly and easily filled using the SAFE-FIL® system.

## Electrical Power

Global electrical supply compatibility, and up to 30 minutes battery back-up make it suitable for use worldwide.

## Low life costs

Minimal service procedures - annual calibration checks, and battery replacement every two years.

## Agent Status and Monitoring

The large capacity reservoir can accommodate more than a whole bottle of Suprane® (desflurane) enabling continuous use for long periods. The illuminated level indicator shows how much agent is available, even in subdued lighting conditions.

### 1 Vaporizer Status Indicators

Visual indication of vaporizer status including battery condition, alarms, and agent level.

### 2 Innovative Filling System

The SAFE-FIL® system prevents accidental release of vapour during refill procedures.

### 3 Vaporizer Control

Easy to use agent control with LED indication when the concentration has been set above 12%.

### 4 Agent Level Indicators

High-intensity LEDs indicate agent level and warn when refilling is required.



1

2

3

4



facebook.com/PenlonGlobal



linkedin.com/company/penlon-ltd



twitter.com/PenlonGlobal

1639

Doc. No. SE0521SL / October 2021. Penlon, InterMed and Sigma EVA are trademarks of Penlon Limited. All other trademarks are the property of their respective owners.



Part of the InterMed Group © Penlon Ltd. 2021 All rights reserved

Penlon Limited  
Abingdon Science Park  
Barton Lane, Abingdon  
OX14 3NB, UK

**General**  
t +44 (0) 1235 547000  
f +44 (0) 1235 547041  
w www.penlon.com

**International Sales**  
t +44 (0) 1235 547001  
f +44 (0) 1235 547021  
e international.sales@penlon.com

**UK Sales**  
t +44 (0) 1235 547036  
f +44 (0) 1235 547023  
e uk.sales@penlon.com

**Technical Support**  
t +44 (0) 1235 547060  
f +44 (0) 1235 547061  
e tech.support@penlon.com

EC REP

OBELIS S.A.  
Bd Général Wahis, 53  
B-1030 Brussels  
Belgium